IN VITRO INTERNATIONAL

Non animal testing

Results available in just hours

**INVITRO IRRITATION
ASSAY SYSTEM**

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Maximal Qualified Score (MQS)	Predicted UN GHS classification**
0 – 12.5	No Category
> 12.5 – 30.0	No Prediction Can be Made*
> 30.0	Category 1

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- | | |
|---------------------------------|-----------------------------|
| ✓ Completely, 100%, animal free | ✓ Results in just 3 minutes |
| ✓ No need for a GLP laboratory | ✓ and no more than 4 hours |
| ✓ No need for cell cultures | ✓ Easy to use |
| ✓ OECD TG, globally accepted | |

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Irritection® Assay System

Toxicology in vitro I risultati di questo studio di convalida sono stati presentati a EURL-ECVAM che lo sta esaminando come un metodo che ha *“potential to be used as a partial replacement to identify chemicals not requiring classification for serious eye damage/eye irritation (as an initial step of a bottom-up testing strategy) as well as chemicals inducing serious eye damage (as an initial step of a top-down testing strategy) according to the UN GHS / EU CLP since its performance seems to be comparable to that of other validated test methods and also because it is a highly cost-effective assay as compared to other methods.”*

DB-ALM Protocol n°157 *“The Ocular Irritection® assay system, represents an optimized protocol of the former Eytex®, and it is a standardized and quantitative in vitro test that predicts the ocular irritation potential of cosmetics, consumer products, pharmaceuticals, and chemical raw materials. The Ocular Irritection® assay system evaluates the ocular hazard effects of test substances based on the premise that eye irritation and corneal opacity after exposure to irritating substances is the result of perturbation or denaturation of corneal proteins.”*

*“These updates have led to use of the Ocular Irritection® assay system over the last 20 years in cosmetics, personal care products, chemical, textile, surfactant, and petrochemical products’ as well as in selected other markets. The assay is currently available, as a kit comprised of reagents and computer software that have been integrated with user instrumentation to provide a standardised and an automated in vitro test. Furthermore, **because of its biochemical nature, Ocular Irritection® kit offers a long shelf-life (years), does not require sterile conditions and is readily available across the world.”***

OECD TG 496: *“The in vitro macromolecular test method is a biochemical in vitro test method that can be used to identify chemicals (substances and mixtures) that have the potential to induce serious eye damage as well as chemicals not requiring classification for eye irritation or serious eye damage.”*

“The in vitro macromolecular test method is recommended as an initial step of a Top-Down testing strategy approach as described within the OECD Guidance Document (GD) 263 (2) to positively identify chemicals inducing serious eye damage, i.e., chemicals to be classified as UN GHS Category 1 (1) without further testing. The in vitro macromolecular test method is also recommended to identify chemicals that do not require classification for eye irritation or serious eye damage as defined by the UN GHS (UN GHS No Category) (1) and may therefore be used as an initial step within a Bottom-Up testing strategy approach (OECD GD 263) (2).”

Corrositex® Assay System

Department of Transportation-DOT-E 10904: Original exemption granted April 28, 1993

Consumer Product Safety Commission (CPSC): Formal Acceptance, NIEHS press release dated 3/21/00

European Centre for the Validation of Alternative Methods (ECVAM): 12/2002

International Air Transportation Association (IATA): Letter of acceptance dated December 17, 1993

National Institute of Environmental Health Sciences (NIEHS): Endorsement dated 6/22/99

Occupational Safety and Health Administration (OSHA): Letter of Interpretation dated March 3, 1994

Formal Acceptance, NIEHS press release dated 3/21/00

Transport Canada - Permit For Equivalent Level of Safety SU 4483: Original approval 8/14/96 - Full Draize Replacement Acceptance 3/5/02

OECD Guideline 435: “In Vitro Membrane Barrier Test Method for Skin Corrosion”, in accordance with **GHS Packing Group**.